

PHILIPS

AVENT

Natural baby bottle

1 Bottle

2 oz/60 ml

First flow teat

0 m



SCF039/17



Easy to combine with breastfeeding Natural latch on

Our Natural bottle with an ultra-soft teat more closely resembles the breast. The wide breast-shaped teat with flexible spiral design and comfort petals allows natural latch-on and makes it easy to combine breast and bottle feeding.

A comfortable and contented feed for your baby

- Flexible spiral design combined with comfort petals
- Anti-colic valve designed to reduce colic and discomfort

Easy to combine breast and bottle feeding

- Natural latch-on due to the wide breast shaped teat
- Ultra soft teat designed to mimic the feel of the breast

Other benefits

- Simple to use and clean, quick and easy assembly
- Teats with different flow rates available

Ideal for newborns and premies

- Smaller 2 oz/60 ml bottle designed for newborn feeding
- Slower-flow teat

Highlights

2-oz/60-ml bottle container

A smaller-size bottle to help ensure the appropriate amount for baby's smaller tummy.

Slower-flow teat

A smaller hole size delivers a more controlled flow rate for babies who drink more slowly.

Natural latch on



The wide breast shaped teat promotes a natural latch-on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Ultra soft teat



The teat has an ultra soft texture, designed to mimic the feel of the breast.

Flexible spiral design



Flexible spiral design, combined with our unique comfort petals to create a flexible teat, allowing for a more natural feed without teat collapse.

Anti-colic valve



Anti-colic valve designed to keep air away from your baby's tummy, to help reduce colic and discomfort.

Choosing the right teat



The Philips Avent Natural range offers different teat softness and increasing flow rates for each of your baby's development stages.

Simple to use and clean



Wide bottle neck makes filling and cleaning easy. Only a few parts for quick and simple assembly.

Specifications

Design

Bottle design: Ergonomic shape, Wide neck

Development stages

Stages: 0-6 months

Ease of use

Bottle use: Easy to assemble, Easy to clean, Easy to hold

Functions

Anti-colic valve: Advanced anti-colic system

Latch on: Easily combine breast and bottle, Natural latch on

Teat: Ultra-soft and flexible teat, unique comfort petals

Material

Bottle: Polypropylene, BPA-free*

Teat: BPA-free*, Silicone

What is included

Baby bottle: 1 pcs



* 0% BPA, following EU regulation 10/2011

* What is colic and how does it affect babies? Colic is caused in part by swallowing air while feeding, which creates discomfort in a baby's digestive system. Symptoms include crying and fussing.