

PHILIPS

AVENT

Natural teat

2 pieces

Medium flow

3m+



SCF043/27



Easy to combine with breastfeeding

Extra soft and flexible teat

Our ultra soft teat with flexible spiral design, more closely resembles the breast. The comfort petals and natural teat shape allows natural latch on and makes it easy to combine breast and bottle feeding.

Easy to combine breast and bottle feeding

- Natural latch on due to the wide breast shaped teat
- Ultra soft teat designed to mimic the feel of the breast

A comfortable and contented feed for your baby

- Flexible spiral design combined with comfort petals

Vents air away from your baby's tummy

- Anti-colic valve designed to reduce colic and discomfort

Other benefits

- Teats with different flow rates available

Highlights

Specifications

Natural latch on



The wide breast shaped teat promotes a natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Flexible spiral design



Flexible spiral design, combined with our unique comfort petals to create a flexible teat, allowing for a more natural feed without teat collapse.

Ultra soft teat



The teat has an ultra soft texture, designed to mimic the feel of the breast.

Anti-colic valve



Anti-colic valve designed to keep air away from your baby's tummy, to help reduce colic and discomfort.

Material

Teat: Silicone, BPA free*

What is included

Ultra soft, medium flow teat: 2 pcs

Development stages

Stages: 3m+

Functions

Latch on: Natural latch on, Easy combine breast and bottle

Teat: Unique comfort petals, Ultra soft and flexible teat

Anti-colic valve

Choosing the right teat



The Philips Avent Natural range offers different teat softness and increasing flow rates for every development stage of your baby.

