

PHILIPS

AVENT

Natural teat

2 pieces

Fast flow

6m+

SCF044/27



Easy to combine with breastfeeding

Natural latch on

Our soft and anti-collapse ribbed teat is designed for growing babies. The comfort petals and natural teat shape allows natural latch-on and makes it easy to combine breast and bottle feeding.

Easy to combine breast and bottle feeding

- Natural latch-on due to the wide breast shaped teat
- Soft and smooth silicone for your baby's changing needs

A comfortable and contented feed for your baby

- Flexible, anti-collapse ribbed teat design

Vents air away from your baby's tummy

- Anti-colic valve designed to reduce colic and discomfort

Other benefits

- Teats with different flow rates available

Highlights

Specifications

Natural latch on



The wide breast shaped teat promotes a natural latch-on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Soft and smooth silicone

Bite resistant and smooth teat designed for your growing baby's changing needs.

Flexible, anti-collapse design



Petals and ridges inside the teat allow flexibility without collapse, for an uninterrupted feed.

Anti-colic valve



Anti-colic valve designed to keep air away from your baby's tummy, to help reduce colic and discomfort.

Choosing the right teat



The Philips Avent Natural range offers different teat softness and increasing flow rates for each of your baby's development stages.

Country of origin

United Kingdom

Material

Teat: Silicone, BPA-free*

What is included

Soft, fast flow teat: 2 pcs

Development stages

Stages: 6m+

Functions

Latch on: Natural latch on, Easily combine breast and bottle

Teat: Unique comfort petals, soft and flexible teat

Anti-colic valve

