

# PHILIPS

## AVENT

### Natural teat

2 pieces

Variable flow

3m+

SCF045/27



## Easy to combine with breastfeeding

### Natural latch on

Our soft and anti-collapse ribbed teat is designed for growing babies. The comfort petals and natural teat shape allows natural latch on and makes it easy to combine breast and bottle feeding.

#### **Easy to combine breast and bottle feeding**

- Natural latch on due to the wide breast shaped teat
- Soft and smooth silicone for your baby's changing needs

#### **A comfortable and contented feed for your baby**

- Flexible, anti-collapse ribbed teat design

#### **Vents air away from your baby's tummy**

- Anti-colic valve designed to reduce colic and discomfort

#### **Other benefits**

- Teats with different flow rates available

# Highlights

## Natural latch on



The wide breast shaped teat promotes a natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

## Soft and smooth silicone

Bite resistant and smooth teat designed for growing baby's changing needs.

## Flexible, anti-collapse design




Petals and ridges inside the teat allow flexibility without collapse, for an uninterrupted feed.

## Anti-colic valve



Anti-colic valve designed to keep air away from your baby's tummy, to help reduce colic and discomfort.

## Choosing the right teat



The Philips Avent Natural range offers different teat softness and increasing flow rates for every development stage of your baby.

# Specifications

## Material

Teat: Silicone, BPA free\*

## What is included

Soft, variable flow teat: 2 pcs

## Development stages

Stages: 3m+

## Functions

**Latch on:** Natural latch on, Easy combine breast and bottle

**Anti-colic valve**

**Teat:** Unique comfort petals, soft and flexible teat

