

# PHILIPS

## AVENT

Natural glass baby bottle

1 Bottle

120ml/4oz

Newborn flow teat

0m+



SCF051/17



## Easy to combine with breastfeeding

### Natural latch on

Our Natural bottle with an ultra soft teat more closely resembles the breast. The wide breast-shaped teat with flexible spiral design and comfort petals allows natural latch on and makes it easy to combine breast and bottle feeding.

#### Easy to combine breast and bottle feeding

- Ultra soft teat designed to mimic the feel of the breast
- Natural latch on due to the wide breast shaped teat
- Flexible spiral design combined with comfort petals

#### Pure glass

- Heat resistance
- Pharmaceutical grade glass

#### Vents air away from your baby's tummy

- Anti-colic valve designed to reduce colic and discomfort

#### Other benefits

- Simple to use and clean, quick and easy assembly
- Teats with different flow rates available
- Compatible with the Philips Avent Natural range
- This bottle is BPA free\*

# Highlights

## Ultra soft teat



The teat has an ultra soft texture, designed to mimic the feel of the breast.

## Natural latch on



The wide breast shaped teat promotes a natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

## Flexible spiral design



Flexible spiral design, combined with our unique comfort petals to create a flexible teat, allowing for a more natural feed without teat collapse.

## Heat resistance



The Natural glass bottle is heat and thermal shock resistant. It can therefore be safely stored in the fridge, warmed, and is also suitable for sterilization.

## Pharmaceutical grade glass

Premium quality borosilicate glass to ensure the highest quality for ultimate purity.

## Anti-colic valve



Anti-colic valve designed to keep air away from your baby's tummy, to help reduce colic and discomfort.

## Simple to use and clean



Wide bottle neck makes filling and cleaning easy. Only a few parts for quick and simple assembly.

## Choosing the right teat



The Philips Avent Natural range offers different teat softness and increasing flow rates for every development stage of your baby.

## Compatibility across the range



Mix and match our breast pump, bottle and cup parts, and create the product that works for you, when you need it!

## This bottle is BPA free\*



The Philips Avent Natural bottle is made of BPA free\* material (polypropylene).

# Specifications

## Material

**Bottle:** Glass, BPA free\*

**Teat:** Silicone, BPA free\*

## What is included

**Baby bottle:** 1 pcs

## Design

**Bottle design:** Ergonomic shape, Wide neck

## Ease of use

**Bottle use:** Easy to assemble, Easy to clean, Easy to hold, Dishwasher & microwave safe

## Development stages

**Stages:** 0-6 months

## Functions

**Latch on:** Natural latch on, Easy combine breast and bottle

**Teat:** Ultra soft and flexible teat, Unique comfort petals

**Anti-colic valve**

