

PHILIPS

AVENT

Classic Pacifier Value
Pack 6-18m, 3 pack

6-18m

Orthodontic & BPA-Free

3-pack

SCF134/30



Orthodontic for maximum comfort

Day & night value pack

Philips Avent orthodontic, collapsible and symmetrical nipples respect the natural development of baby's palate, teeth and gums. All Philips Avent pacifiers are made of silicone and are taste and odor-free. Colors are subject to change.

Orthodontic nipple

- Orthodontic, symmetrical collapsible nipple

Hygiene

- Snap-on protective cap
- Can be sterilized
- Dishwasher Safe

Safety

- Security ring handle

Odorless and taste-free

- User-friendly silicone nipples

Unique security ring handle

- Glow in the dark handle

Highlights

Orthodontic nipple



Philips Avent flat, drop-shaped symmetrical nipples respect the natural development of your baby's palate, teeth and gums, even if the pacifier ends upside-down in the mouth.

Security ring handle



For easy removal of the Philips Avent pacifier at any time

Snap-on protective cap



To keep sterilized nipples hygienic

User-friendly silicone nipples



The Philips Avent silicone nipple is taste and odor-free so it's more likely to be accepted by your baby. The silicone is smooth, transparent, easy to clean and it doesn't get sticky. The nipple is strong, long-lasting and won't become mis-shapen or discolored over time.

Can be sterilized



Can be sterilized

Dishwasher Safe



Dishwasher Safe

Glow in the dark handle



A night pacifier with a unique security handle that glows in the dark.

Specifications

Country of origin

England

What is included

Silicone Pacifier: 3 pcs

Snap-on protective cap: 3 pcs

Features

Can be sterilized

BPA-Free

Helps to comfort your baby

Security ring handle

Development stages

Stages: 6-18 months

