

# PHILIPS

## AVENT

Milk powder  
dispenser

3 doses



SCF135/06



## Ideal for travel

BPA Free

This unit carries 3 pre-measured portions of milk powder in separate compartments. When you're ready to feed, just pour the powder into the feeding bottle of pre-boiled cooled water. Remove inner sections to use as a bowl or container.

**Benefits and features:**

- Holds enough powdered formula for three 260 ml/ 9 oz feeds
- Entire dispenser sterilizable, microwavable & dishwasher safe
- Inner part can be removed

# Highlights

## Converts into handy snack cup



Inner part can be removed to convert into a handy snack cup

## Holds enough formula for three 260 ml/ 9 oz feeds



The Philips Avent milk powder dispenser holds 3 pre-measured doses of milk powder - ideal for travel

## Sterilizable, microwavable, dishwasher safe

All parts can be sterilized, microwaved and are dishwasher safe for quick and easy cleaning

# Specifications

## What is included

Milk powder dispenser: 1 pcs

## Country of origin

Malaysia

## Material

BPA free\*

## Development stages

Stages: 0-6 months, 0 - 6 months

