

PHILIPS

AVENT

Milk powder
dispenser

3 doses



SCF135/06



Ideal for travel

BPA-Free

This handy little unit carries 3 pre-measured portions of milk powder in separate compartments. When you're ready to feed, just pour the powder into the bottle of pre-boiled cooled water. Remove inner sections to use as a bowl or container.

Benefits and features:

- Holds enough powdered formula for three 260 ml/ 9 oz feeds
- Entire dispenser sterilizable, microwavable & dishwasher safe
- Inner part can be removed

