

PHILIPS

AVENT

Milk powder
dispenser

3 doses

SCF135/18



Ideal for travel

BPA Free

This unit carries 3 pre-measured portions of milk powder in separate compartments. When you're ready to feed, just pour the powder into the feeding bottle of pre-boiled cooled water. Remove inner sections to use as a bowl or container.

Benefits and features:

- Inner part can be removed
- Entire dispenser sterilizable, microwavable & dishwasher safe
- Holds enough powdered formula for three 260 ml/ 9 oz feeds

Highlights

Converts into handy snack cup



Inner part can be removed to convert into a handy snack cup

Holds enough formula for three 260 ml/ 9 oz feeds



The Philips Avent milk powder dispenser holds 3 pre-measured doses of milk powder - ideal for travel

Sterilizable, microwavable, dishwasher safe

All parts can be sterilized, microwaved and are dishwasher safe for quick and easy cleaning

Specifications

Country of origin

Malaysia

Development stages

Stages: 0 - 6 months, 0-6 months

Material

BPA free*

What is included

Milk powder dispenser: 1 pcs

