



Milk powder dispenser

- 3 doses

SCF135/18

## Ideal for travel

BPA Free

This unit carries 3 pre-measured portions of milk powder in separate compartments. When you're ready to feed, just pour the powder into the feeding bottle of pre-boiled cooled water. Remove inner sections to use as a bowl or container.

### Benefits

#### Benefits and features:

- Inner part can be removed
- Holds enough powdered formula for three 260 ml/ 9 oz feeds
- Entire dispenser sterilizable, microwavable & dishwasher safe

### Features

#### Converts into handy snack cup



Inner part can be removed to convert into a handy snack cup

**PHILIPS**

Holds enough formula for three 260 ml/ 9 oz feeds



The Philips Avent milk powder dispenser holds 3 pre-measured doses of milk powder - ideal for travel

Sterilizable, microwavable, dishwasher safe

All parts can be sterilized, microwaved and are dishwasher safe for quick and easy cleaning

## Specifications

Please note that this is a pre-sales leaflet. The contents of this leaflet reflect the best of our knowledge per date and country mentioned above. The contents of this leaflet are subject to change without notice. Philips does not accept any liability as to the contents of this leaflet.

### Country of origin

Malaysia Yes

### Material

BPA free\* Yes

### Development stages

Stages 0 - 6 months  
0-6 months

### What is included

Milk powder dispenser 1 pcs



data subject to change  
2018, October 21  
Version: 8.0.1  
UPC: 00075020046185

© 2018 Koninklijke Philips N.V.  
All rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

[www.philips.com](http://www.philips.com)