

**PHILIPS**

**AVENT**

Nipple travel pack

SCF139/02



## Keeps Avent teats sterile for up to 24 hours

Ideal nipple for travel

Philips Avent teat travel pack SCF139/02 is extremely useful to keep teats sterile until you need them to feed your baby.

**Benefits and features:**

- Ideal for travel

# Specifications

## Country of origin

England

## What is included

Dome cap: 2 pcs

Screw ring: 2 pcs

Sealing base: 2 pcs

## Development stages

Stages: 0-6 months, 0 - 6 months, 6 - 12 months

