

# PHILIPS

## AVENT

Sealing discs for  
feeding bottle

Bottle accessories

SCF143/06



## Easy milk and food storage

Ideal for storing and freezing breast milk or food

The Philips AVENT sealing discs convert all Philips AVENT feeding bottles to breast milk and food storage containers. They are compatible with our Natural and Classic bottles and BPA free.

### **Easily store your milk and food**

- 6 sealing discs to store breast milk or food
- Ideal for freezing breast milk
- Write date on the disc for storage

### **0% BPA (Bisphenol A Free)**

- Philips AVENT sealing discs are BPA free

# Highlights

## 6 sealing discs


Easily store your milk and food in the Philips Avent bottle. Close the bottle with the sealing disc.

## Ideal for freezing breast milk



The sealing disc allow you to store or freeze your precious breast milk, without leaking. Please keep in mind that the sealing disc has to be used with the Philips Avent adapter ring in order to prevent leaking.

## Write date on the disc



You can write the date of expressing or freezing on the sealing disc. Easily organize your breast milk or baby food.

## BPA free disc



Philips Avent sealing discs are made from polypropylene (PP) & thermoplastic elastomer (TPE), two BPA free materials.

# Specifications

## Country of origin

England

## Material

Thermoplastic Elastomer (TPE)  
BPA free\*

## Development stages

Stages: 0-6 months

## Functions

**Ease of use:** Write date on disc, Non leakage

**Easy storage:** Store/freeze milk and food

## What is included

Sealing disc for milk storage: 6 pcs

