

PHILIPS

AVENT

Sealing discs for
feeding bottle

Bottle accessories

SCF143/06



Easy milk and food storage

Ideal for storing and freezing breast milk or food

The Philips AVENT sealing discs convert all Philips AVENT feeding bottles into breast milk and food storage containers. They are compatible with our Natural and Classic bottles and are BPA-free.

Easily store your milk and food

- 6 sealing discs to store breast milk or food
- Ideal for freezing breast milk
- Write date on the disc for storage

BPA Free (0% Bisphenol A)

- Philips AVENT sealing discs are BPA free

Highlights

6 sealing discs

Easily store your milk and food in the Philips Avent bottle. Close the bottle with the sealing disc.

Ideal for freezing breast milk



The sealing disc allows you to store or freeze your precious breast milk, without leaking. Please keep in mind that the sealing disc has to be used with the Philips Avent adapter ring in order to prevent leaking.

Write date on the disc



You can write the date of expressing or freezing on the sealing disc. Easily organise your breast milk or baby food.

BPA free disc



Philips Avent sealing discs are made from polypropylene (PP) and thermoplastic elastomer (TPE), two BPA-free materials.

Specifications

Country of origin

England

Material

Thermoplastic Elastomer (TPE)
BPA free*

Development stages

Stages: 0-6 months

Functions

Ease of use: Write date on disc, Non-leakage

Easy storage: Store/freeze milk and food

What is included

Sealing disc for milk storage: 6 pcs

