

PHILIPS

AVENT

Mini soother

Perfect fit for your newborn

0-2m

Orthodontic and BPA-Free

2-pack



SCF151/03



Extra small and lightweight for your little one A comfortable fit for your newborn

The Philips AVENT Mini Soother soothes tiny babies up to 2 months old. Our small, lightweight shield comfortably fits your newborn without touching their nose. Our orthodontic collapsible teat respects your little one's oral development.

For tiny babies

- Extra-small and lightweight for tiny babies

Comfortable soft silicone teat

- 9 out of 10 babies accept our soothers*
- Designed for natural oral development

Quality assured

- Made at our award-winning site in the UK

Safety and Hygiene

- Security handle for easy removal
- Snap-on cap helps keep your baby's soother clean

Easy to clean

- Easy to sterilise for extra hygiene

Highlights

Extra small and lightweight



Mini Soother's shield is extra-small and lightweight for tiny babies up to 2 months old.

Teat acceptance



Babies know what they like! We asked mums how their little ones respond to Philips AVENT teats and 9 out of 10 babies accept our soothers.*

Orthodontic teat



Our collapsible silicone nipple has a symmetrical shape that respects your baby's palate, teeth and gums as he or she grows.

Made in the UK



You can be sure that your little one's comfort is in good hands. This soother was made at our award-winning site in the UK.*

Security handle



Our security handle lets you easily remove your baby's soother at any time. Even little hands can grab it!

Snap-on cap for hygiene



When the soother isn't in use, simply snap on the cap before storing to keep the teat safe and clean.

Easy to sterilise



Keeping your little one's Soothies and soothers clean is easy. Simply put them in your steriliser or submerge in boiling water.

Specifications

What is included

Mini soother: 2 pcs

Safety

Safety ring handle
BPA-free

Hygiene

Easy to clean
Can be sterilised
Dishwasher safe

