

Philips Avent
Nipple Protector

Small

SCF156/00



Helping you to breastfeed longer

Breast care, protects sore nipples

Philips Avent nipple protector SCF156/00 made of ultra-fine, soft, odourless, taste-free silicone that protect sore or cracked nipples during breastfeeding.

Avert sore and cracked nipples when breast feeding

- Protect sore nipples during breastfeeding

Helps you to breast feed longer

- Easy latch on for your baby
- Made with odorless, taste free, ultra-fine silicone
- Shaped to allow skin contact with baby

PHILIPS
AVENT

Nipple Protector
Small

SCF156/00

Specifications

Country of origin

- England

What is included

- Small Nipple Protector: 2 pcs

Weight and dimensions

- Diameter: 15 mm

Development stages

- Stages: 0 - 6 months

Design

- Allows baby to feel your skin
- Allows baby to smell your skin

Ease of use

- Heal when still breast feeding

Functions

- Easy latch on: Ultra fine silicone

Material

- Nipple protector: BPA Free*, Silicone

Highlights

Shaped to allow skin contact with baby

Your baby can still feel and smell your skin and continue to stimulate your milk supply while suckling, and will return easily to the breast once your nipples are healed.

Protect sore nipples during breastfeeding

Philips Avent Nipple Protectors are only designed for use when you have sore or cracked nipples and should be used with health professional advice.



Issue date 2018-07-03

Version: 14.3.1

12 NC: 0000 000 00000

© 2018 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com