

**PHILIPS**

**AVENT**

Nipple Protector

Medium(21mm)

2 pcs

SCF156/01



## Helping you to breastfeed longer

### Breast care, protects sore nipples

Philips Avent nipple protector SCF156/00 made of ultra-fine, soft, odourless, taste-free silicone that protect sore or cracked nipples during breastfeeding.

#### **Avert sore and cracked nipples when breast feeding**

- Protect sore nipples during breastfeeding

#### **Helps you to breast feed longer**

- Easy latch on for your baby
- Made with odorless, taste free, ultra-fine silicone
- Shaped to allow skin contact with baby

# Highlights

## Protect sore nipples during breastfeeding

Philips Avent Nipple Protectors are only designed for use when you have sore or cracked nipples and should be used with health professional advice.

## Easy latch on for your baby

Your baby can easily latch on through the shield and create a seal.

## Made of ultra-fine silicone

The Philips Avent Nipple Protectors are made of an odorless, tasteless ultra-fine silicone.

## Shaped to allow skin contact with baby

Your baby can still feel and smell your skin and continue to stimulate your milk supply while suckling, and will return easily to the breast once your nipples are healed.

# Specifications

## Country of origin

United Kingdom

## What is included

Standard Nipple Protector: 2 pcs

## Weight and dimensions

Diameter: 21 mm

## Development stages

Stages: 0-6 months

## Design

Allows baby to feel your skin  
Allows baby to smell your skin

## Ease of use

Heal when still breast feeding

## Functions

Easy latch on: Ultra fine silicone

## Material

Nipple protector: BPA Free\*, Silicone

