

PHILIPS

AVENT

Comfort breast shell
set

2 pcs



SCF157



Comfort and protection

Ultra comfortable breast shells

The ultrasoft Philips Avent breast shells SCF157/02 are worn inside your bra to protect your nipples from chafing and to collect excess breast milk.

Maximum comfort

- Collect excess breast milk
- Help ease engorgement

Maximum protection

- Protect sore nipples
- The breast shells can be used when breast feeding or pumping



asimpleswitch.com

Highlights

Collect excess breast milk

Philips Avent breast milk collection shells (no holes)- collect excess breast milk when feeding or using a breast pump.

Protect sore nipples

Ventilated shells- protect sore or cracked nipples to help them heal more quickly. Their gentle pressure helps relieve engorgement. The holes allow air to circulate.

Specifications

Material

Backing cushions: Silicone

What is included

Ventilated breast shells: 2 pcs

Non-ventilated breast milk saver shells: 2 pcs

Ultrasoft backing cushions: 2 pcs

Development stages

Stages: 0 - 6 months

Protection against

Cracked nipples

Sore nipples



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

* This section contains consumer opinions on the product. Philips dissociates itself from the contents entered by consumers in this section and consequently any technical information and / or advice on the use of the product included therein are not intended as official Philips information.