

Classic soother

- For essential comfort
- 0–6m
- Orthodontic and BPA-Free
- 2-pack

SCF169/26

Designed for your baby's everyday comfort needs

Soother with bright, colourful designs

Help meet your little one's essential soothing needs around the clock with the Philips Avent Classic soother. Available in a range of colours, our orthodontic collapsible teat respects your baby's natural oral development.

Benefits

Comfortable soft silicone teat

- 9 out of 10 babies accept our soothers*
- Designed for natural oral development

Quality assured

- Made at our award-winning site in the UK

Safety and Hygiene

- Security handle for easy removal
- Snap-on cap helps keep your baby's soother clean

Easy to clean

- Easy to sterilise for extra hygiene

PHILIPS

Features

Teat acceptance



Babies know what they like! We asked mums how their little ones respond to Philips AVENT teats and 9 out of 10 babies accept our soothers.*

Orthodontic teat




Our collapsible silicone nipple has a symmetrical shape that respects your baby's palate, teeth and gums as he or she grows.

Made in the UK




You can be sure that your little one's comfort is in good hands. This soother was made at our award-winning site in the UK.*

Security handle



Our security handle lets you easily remove your baby's soother at any time. Even little hands can grab it!

Snap-on cap for hygiene



When the soother isn't in use, simply snap on the cap before storing to keep the teat safe and clean.

Easy to sterilise



Keeping your little one's Soothies and soothers clean is easy. Simply put them in your steriliser or submerge in boiling water.

Specifications

Please note that this is a pre-sales leaflet. The contents of this leaflet reflect the best of our knowledge per date and country mentioned above. The contents of this leaflet are subject to change without notice. Philips does not accept any liability as to the contents of this leaflet.

What is included

Classic soother 2 pcs

Safety

Safety ring handle Yes
BPA-free Yes

Hygiene

Easy to clean Yes
Can be sterilised Yes
Dishwasher safe Yes

* No. 1 global soother brand

* For hygiene reasons, replace soothers after 4 weeks of use

* Our range supports mums and babies at every development stage

* Tested online with 100 mums, UK 2012

* 2014 Manufacturer of the Year



data subject to change
2019, January 17

Version: 3.0.1
EAN: 08710103811091

© 2019 Koninklijke Philips N.V.
All rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

www.philips.com