

PHILIPS

AVENT

Translucent
Soothers

0-3m

BPA-Free

SCF170/20



Orthodontic BPA-Free

Stylish design, crystal clear look.

Philips AVENT orthodontic, collapsible and symmetrical teats SCF170/20 respect the natural development of baby's palate, teeth and gums. All AVENT soothers are made of silicone and are taste- and odour-free. Colours are subject to change.

Orthodontic teat

- Orthodontic, symmetrical collapsible teat

Odourless and taste-free

- User-friendly silicone teats

Hygiene

- Snap-on protective cap

Safety

- Security ring handle

Highlights

Orthodontic teat



Philips AVENT flat, drop-shaped symmetrical teats respect the natural development of your baby's palate, teeth and gums, even if the soother ends up upside down in the mouth.

User-friendly silicone teats



The Philips AVENT silicone teat is taste- and odour-free so it's more likely to be accepted by your baby. The silicone is smooth, transparent


and easy to clean, and it doesn't get sticky. The teat is strong and long-lasting, and won't become misshapen or discoloured over time.

Snap-on protective cap



To keep sterilised teats hygienic

Security ring handle



For easy removal of the Philips AVENT soother at any time

Specifications

Country of origin

England

What is included

Silicone Soother: 2 pcs

Snap-on protective cap: 2 pcs

Development stages

Stages: 0–6 months

