

PHILIPS

AVENT

Night-time soother

With glow-in-the-dark handle

0–6m

Orthodontic and BPA-Free

2-pack



SCF176/28



Convenient soothing at bedtime

Unique glow-in-the-dark handle

Comfort your baby at bedtime with the Philips Avent Classic Night-Time soother: our glow-in-the-dark handle makes it easier to find at night. As they sleep, our orthodontic collapsible teat respects your baby's natural oral development.

Glow-in-the-dark handle*

- Find glowing handle more easily at bedtime

Comfortable silicone teat

- 9 out of 10 babies accept our soothers*
- Designed for natural oral development

Quality assured

- Made at our award-winning site in the UK

Safety and Hygiene

- Security handle for easy removal
- Snap-on cap helps keep your baby's soother clean

Easy to clean

- Easy to sterilise for extra hygiene

Highlights

Glow-in-the-dark handle



We know that getting your little one back to sleep is important. Our unique glow-in-the-dark handle helps you to find this soother without having to switch on the lights.*

Teat acceptance



Babies know what they like! We asked mums how their little ones respond to Philips AVENT teats and 9 out of 10 babies accept our soothers.*

Orthodontic teat



Our collapsible silicone nipple has a symmetrical shape that respects your baby's palate, teeth and gums as he or she grows.

Made in the UK



You can be sure that your little one's comfort is in good hands. This soother was made at our award-winning site in the UK.*

Security handle




Our security handle lets you easily remove your baby's soother at any time. Even little hands can grab it!

Snap-on cap for hygiene



When the soother isn't in use, simply snap on the cap before storing to keep the teat safe and clean.

Easy to sterilise



Keeping your little one's Soothies and soothers clean is easy. Simply put them in your steriliser or submerge in boiling water.

Specifications

What is included

Night-time classic soother: 2 pcs

Safety

Safety ring handle
BPA-free

Hygiene

Easy to clean
Can be sterilised
Dishwasher safe

