

PHILIPS

AVENT

Contemporary
Freeflow Pacifier

Soothers

0-6m

BPA-Free



SCF180/23



Extra airflow for baby's skin

Freeflow shields and modern contemporary designs.

Philips AVENT orthodontic, collapsible teats SCF180/23 respect the natural development of baby's palate, teeth and gums. All our soothers are made from silicone and are taste and odor-free. Colors are subject to change.

Orthodontic nipple

- Orthodontic, symmetrical collapsible nipple

Odorless and taste-free

- User-friendly silicone nipples

Hygiene

- Snap-on protective cap
- Can be sterilized
- Dishwasher Safe

Safety

- Security ring handle



asimpleswitch.com

Highlights

Orthodontic nipple




Philips AVENT flat, drop-shaped symmetrical nipples respect the natural development of your baby's palate, teeth and gums, even if the pacifier ends upside-down in the mouth.

User-friendly silicone nipples



The Philips AVENT silicone teat is taste and odor-free so it's more likely to be accepted by your baby. The silicone is smooth, transparent, easy to clean and it doesn't get sticky. The nipple is strong, long-lasting and won't become mis-shaped or discolored over time.

Security ring handle



For easy removal of the Philips AVENT pacifier at any time

Can be sterilized



Can be sterilized

Dishwasher Safe



Dishwasher Safe

Snap-on protective cap



To keep sterilized nipples hygienic

Specifications

Country of origin

England

What is included

Silicone Pacifier: 2 pcs
Snap-on hygienic cap: 2 pcs

Features

Can be sterilized
BPA-Free
Helps to comfort your baby
Hygienic cap

Development stages

Stages: 0 - 6 months

