

# PHILIPS

## AVENT

Contemporary  
Freeflow Pacifier

Soothers

6-18m

BPA-Free



SCF180/24



## Extra airflow for baby's skin

Freeflow shields and modern contemporary designs.

Philips AVENT orthodontic, collapsible and symmetrical nipples SCF180/24 respect the natural development of baby's palate, teeth and gums. All AVENT pacifiers are made of silicone and are taste and odor-free. Colors are subject to change.

### Orthodontic nipple

- Orthodontic, symmetrical collapsible nipple

### Odorless and taste-free


- User-friendly silicone nipples

### Hygiene

- Snap on hygienic cap
- Can be sterilized
- Dishwasher Safe

### Safety

- Security ring handle



# Highlights

## Orthodontic nipple




Philips AVENT flat, drop-shaped symmetrical nipples respect the natural development of your baby's palate, teeth and gums, even if the pacifier ends upside-down in the mouth.

## User-friendly silicone nipples



The Philips AVENT silicone teat is taste and odor-free so it's more likely to be accepted by your baby. The silicone is smooth, transparent, easy to clean and it doesn't get sticky. The nipple is strong, long-lasting and won't become mis-shaped or discolored over time.

## Security ring handle



For easy removal of the Philips AVENT pacifier at any time

## Snap on hygienic cap



To keep sterilized nipples hygienic.

## Can be sterilized



Can be sterilized

## Dishwasher Safe



Dishwasher Safe

# Specifications

## Country of origin

England

## What is included

Silicone Pacifier: 2 pcs

Snap-on hygienic cap: 2 pcs

## Features

Can be sterilized

BPA-Free

Helps to comfort your baby

Hygienic cap

## Development stages

Stages: 6-18 months

