

# PHILIPS

## AVENT

### Classic soother

For essential comfort

0–6m

Orthodontic and BPA-Free

2 pack

SCF182/12



## Designed for your baby's everyday comfort needs

### Soother with colourful animal designs

Help meet your little one's essential soothing needs around the clock with the Philips Avent Classic soother. Available in a range of colours, our orthodontic collapsible teat respects your baby's natural oral development.

#### Comfortable soft silicone teat

- 9 out of 10 babies accept our soothers\*
- Designed for natural oral development

#### Quality assured

- Made at our award-winning site in the UK

#### Safety and Hygiene

- Security handle for easy removal
- Snap-on cap helps keep your baby's soother clean

#### Easy to clean

- Easy to sterilise for extra hygiene

# Highlights

## Teat acceptance



Babies know what they like! We asked mums how their little ones respond to Philips Avent teats and 9 out of 10 babies accept our soothers.\*

## Orthodontic teat



Our collapsible silicone nipple has a symmetrical shape that respects your baby's palate, teeth and gums as he or she grows.

## Made in the UK



You can be sure that your little one's comfort is in good hands. This soother was made at our award-winning site in the UK.\*

## Security handle




Our security handle lets you easily remove your baby's soother at any time. Even little hands can grab it!

## Snap-on cap for hygiene



When the soother isn't in use, simply snap on the cap before storing to keep the teat safe and clean.

## Easy to sterilise



Keeping your little one's Soothies and soothers clean is easy. Simply put them in your steriliser or submerge in boiling water.

# Specifications

## What is included

Classic soother: 2 pcs

## Safety

Safety ring handle  
BPA-free

## Hygiene

Easy to clean  
Can be sterilised  
Dishwasher safe

