

# PHILIPS

## AVENT

Classic Pacifier 0-6m, 2 pack

For essential comfort

0-6m

Orthodontic & BPA-Free

2-pack

SCF182/13



## Designed for your baby's everyday comfort needs Pacifier with colorful animal designs

Help meet your little one's essential soothing needs around the clock with the Philips Avent Classic pacifier. Available in a range of colors, our orthodontic collapsible nipple respects your baby's natural oral development.

### Comfortable soft silicone nipple

- 9 out of 10 babies accept our pacifiers\*
- Designed for natural oral development

### Quality assured

- Made at our award-winning site in the UK

### Safety & Hygiene

- Security handle for easy removal
- Snap-on cap helps keep your baby's pacifier clean

### Easy to clean

- Easy to sterilize for extra hygiene

# Highlights

## Nipple acceptance



Babies know what they like! We asked moms how their little ones respond to Philips Avent nipples and 9 out of 10 babies accept our pacifiers.\*

## Orthodontic nipple



Our collapsible silicone nipple has a symmetrical shape that respects your baby's palate, teeth and gums as he or she grows.

## Made in the UK




You can be sure that your little one's comfort is in good hands. This pacifier was made at our award-winning site in the UK.\*

## Security handle



Our security handle lets you easily remove your baby's pacifier at any time. Even little hands can grab it!

## Snap-on cap for hygiene



When the pacifier isn't in use, simply snap on the cap before storing to keep the nipple safe and clean.

## Easy to sterilize



Keeping your little one's Soothies and pacifiers clean is easy. Simply put them in your sterilizer or submerge in boiling water.

# Specifications

## What is included

Classic pacifier: 2 pcs

## Safety

Safety ring handle  
BPA-free

## Hygiene

Easy to clean  
Can be sterilized  
Dishwasher safe

