

PHILIPS

AVENT

Classic pacifier 6-18m, 2 pack

For essential comfort

6-18m

Orthodontic & BPA-Free

2-pack



SCF182/15



Designed for your baby's everyday comfort needs

Pacifier for maximum comfort. Fun range

Help meet your little one's essential soothing needs around the clock with the Philips Avent Classic pacifier. Available in a range of colors, our orthodontic collapsible nipple respects your baby's natural oral development.

Comfortable soft silicone nipple

- Designed for natural oral development

Quality assured

- Made at our award-winning site in the UK

Safety & Hygiene

- Security handle for easy removal
- Snap-on cap helps keep your baby's pacifier clean

Easy to clean

- Easy to sterilize for extra hygiene

Highlights

Orthodontic nipple



Our collapsible silicone nipple has a symmetrical shape that respects your baby's palate, teeth and gums as he or she grows.

Made in the UK



You can be sure that your little one's comfort is in good hands. This pacifier was made at our award-winning site in the UK.*

Security handle



Our security handle lets you easily remove your baby's pacifier at any time. Even little hands can grab it!

Snap-on cap for hygiene



When the pacifier isn't in use, simply snap on the cap before storing to keep the nipple safe and clean.

Easy to sterilize



Keeping your little one's Soothies and pacifiers clean is easy. Simply put them in your sterilizer or submerge in boiling water.

Specifications

What is included

Classic pacifier: 2 pcs

Safety

Safety ring handle
BPA-free

Hygiene

Easy to clean
Can be sterilized
Dishwasher safe

