

# PHILIPS

## AVENT

Advanced  
orthodontic pacifiers

6-18m

SCF184/14



## Designed to help healthy oral development

Developed with leading orthodontist Dr. Hagemann

Philips Avent advanced orthodontic teat SCF184/14 is designed to help healthy oral development. All our pacifiers are made from silicone and are taste and odor-free. Colors are subject to change.

### Orthodontic nipple

- Unique "wings"
- Shaped nipple

### Hygiene

- Snap on hygienic cap
- Can be sterilized
- Dishwasher Safe

### Odorless and taste-free

- User-friendly silicone nipples

### IDHF approved

- This product is approved by IDHF

# Highlights

## Unique "wings"



Advanced orthodontic teat with unique wings that minimizes pressure on gums and developing teeth. The wings make the teat wider so that the pressure that is caused by baby's suckling is distributed more evenly resulting in less pressure per tooth.

## Shaped nipple



The nipple is shaped to let baby's tongue remain in natural position.

## Snap on hygienic cap



To keep sterilized nipples hygienic.

## User-friendly silicone nipples

The Philips Avent silicone nipple is taste and odor-free so it's more likely to be accepted by your baby. The silicone is smooth, transparent, easy to clean and it doesn't get sticky. The nipple is strong, long-lasting and won't become mis-shapen or discolored over time.

## Can be sterilized



Can be sterilized

## Dishwasher Safe



Dishwasher Safe

# Specifications

## Country of origin

England

## What is included

Silicone pacifier: 1 pcs

Snap-on hygienic cap: 1 pcs

## Features

Can be sterilized

BPA-Free

Helps to comfort your baby

Hygienic cap

## Development stages

Stages: 6-18 months

