

PHILIPS

AVENT

Soothie pacifier

3m+

Orthodontic & BPA-Free

1 pcs

SCF193/00



Pacifier distributed in Hospitals. BPA-Free

Used by medical professionals to calm newborns

Over 2,000 hospitals in the US give newborns the Philips Avent Soothie SCF190/01 pacifier, which are sized and shaped for baby's developing mouth. Its one-piece construction adheres to the American Academy of Pediatric guidelines.

Hospital grade silicone

- Highly durable

Orthodontic nipple

- Promotes natural development of teeth and gums

Hygienic

- Can be sterilized in a Philips Avent Sterilizer
- Dishwasher safe

BPA Free

- 0% BPA Product

Highlights

0% BPA



This product contains 0% BPA

Can be sterilized

Can be sterilized in a Philips Avent Sterilizer

Dishwasher safe



Dishwasher safe

Highly durable

Highly durable

Promotes natural development

Specially designed to ensure natural development of baby teeth and gums

Specifications

Country of origin

USA

What is included

Silicone Pacifier: 1 pcs

Features

Can be sterilized

BPA-Free

Helps to comfort your baby

Medical grade silicone

One-piece construction

Development stages

Stage: 3 months +

