

PHILIPS

AVENT

Teether Elephant

3m+

Front, middle and back teeth



SCF199/00



Helps soothe baby's gums during teething Reaches front, middle and back teeth

The Philips Avent teether is designed to help soothe teething pain. The teether can be cooled in the fridge, which makes it even more soothing for your little one's sore gums during teething

Choice of various textures

- Soothes and massages sore gums

Easy to hold

- Lightweight design makes it easy for babies to grasp

Easy to clean

- Simply rinse under water or sterilize

Suitable for cooling

- Cool in fridge for extra comfort

BPA Free

- Philips Avent teethers are completely BPA & Phthalate Free



asimpleswitch.com

Highlights

Soothes and massages sore gums

Soothes and massages sore gums

Lightweight design

Lightweight design makes it easy for babies to grasp

Cool in fridge

Cool in fridge for extra comfort

Rinse under water or sterilize



Simply rinse under water or sterilize

BPA & Phthalate Free



BPA free - following EU Directive 2011/8/EU.

Specifications

Country of origin

Malaysia

Features

Helps soothe teething pain

Hygienic

Can be sterilized

BPA-Free

Can be cooled

Development stages

Stages: 3 months +, For front, middle and back

Material

Textures and materials: 2



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.