

PHILIPS

AVENT

PP Adapter ring for
feeding bottles

Classic

SCF200/00



Adapter ring for Classic bottle

Avent feeding bottle

The adapter ring should always be used with the Philips Avent Classic feeding bottle as part of the Classic two-piece anti-colic system.

Benefits and features:

- 0% BPA Product

Two piece anti-colic system

- Allows air into the bottle instead of your baby's tummy

Highlights

Unique anti-colic valve



The two piece Philips Avent anti-colic system consists of the Avent Nipple and adapter ring. As your baby feeds, the unique skirt on the

Avent Nipple flexes to allow air into the baby bottle instead of your baby's tummy. Like breastfeeding, your baby controls the milk flow.

0% BPA



This product contains 0% BPA

Specifications

Country of origin

England

What is included

PP Adapter Ring: 3 pcs

Development stages

Stages: 0-6 months

Material

BPA free*

