

PHILIPS

AVENT

ultra soft pacifier

Ultra soft and flexible

6-18m

Orthodontic & BPA-Free

2-pack



SCF212/23



The softest soother for your baby's delicate skin*

Flexible shield for a comfortable fit

Take care of delicate skin with the Philips Avent ultra soft pacifier. Our super soft, flexible shield follows the shape of your baby's cheeks, leaving fewer marks and causing less skin irritation for more comfortable soothing.

FlexiFit shield technology

- Soft, flexible shield for fewer skin marks & less irritation
- Rounded shield for extra comfort every day
- Extra air holes let your baby's skin breathe

Comfortable soft & silky nipple

- 98% of babies accept ultra soft textured silicon nipples*
- Textured, silky nipple for soothing comfort
- Designed for natural oral development

Sterilisation & carrying case

- Sterilizing in the microwave is safe and easy
- Simple steps for reliable sterilization in 3 minutes

Quality assured

- Made at our award-winning site in the UK



asimpleswitch.com

Highlights

Flexible, soft shield




Delicate skin needs extra care. Our FlexiFit technology lets this shield follow the natural curves of your baby's face for a comfortable fit. Your little one will experience fewer skin marks and less irritation.

Rounded shield



Our rounded shield minimizes pressure for comfortable soothing that's gentle on your little one's cheeks.

Extra airflow



Skin needs to breathe, especially your little one's. Our shield has 6 air holes for extra airflow, designed to reduce skin irritation.

98% accept silicon nipples



Babies know what they like! We asked moms how their little ones respond to our textured silicon nipples and 98% said that their baby accepts the Philips Avent ultra soft pacifiers.

Silky textured nipple



This pacifier's silky, textured silicone nipple soothes and comforts your little one.

Orthodontic nipple



Our collapsible silicone nipple has a symmetrical shape that respects your baby's palate, teeth and gums as he or she grows.

Sterilizer / carrying case



With our carry case you can sterilize your baby's pacifier in the microwave and store hygienically when you're on the move.

Case makes sterilizing easy



Place well-rinsed pacifiers into the case, add 25ml/9fl.oz. of water and microwave at 750-1000W for 3 minutes. Leave to cool for 5 minutes, pour water away and they're ready to go.

Made in the UK



You can be sure that your little one's comfort is in good hands. This pacifier was made at our award-winning site in the UK.*



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

Specifications

What is included

ultra soft pacifier: 2 pcs

Accessories included

Sterilizer/carrying case

Safety

Safety ring handle
BPA-free

Hygiene

Easy to clean
Can be sterilized
Dishwasher safe



- * Developed with healthcare professionals and moms
- * 2016-2017 US consumer tests show an average of 98% teat acceptance of the textured Philips Avent teat used in our ultra air and ultra soft soothers
- * 2014 Manufacturer of the Year
- * For hygiene reasons, replace pacifiers after 4 weeks of use
- * No 1 global pacifier brand
- * Our range supports moms and babies at every development stage
- * 85% of moms surveyed perceive this pacifier feels softer than eight comparable models of other leading brands, independent research, USA, Feb 2017