

#### Ultra Soft Pacifier 6-18m, 2 pack

- Ultra soft and flexible
- 6-18m
- Orthodontic & BPA-Free
- 2-pack

SCF213/22

# The softest soother for your baby's delicate skin\*

Flexible shield for a comfortable fit

Take care of delicate skin with the Philips Avent ultra soft pacifier. Our super soft, flexible shield follows the shape of your baby's cheeks, leaving fewer marks and causing less skin irritation for more comfortable soothing.

#### Benefits

##### FlexiFit shield technology

- Soft, flexible shield for fewer skin marks & less irritation
- Rounded shield for extra comfort every day
- Extra air holes let your baby's skin breathe

##### Comfortable soft & silky nipple

- 98% of babies accept ultra soft textured silicone nipples\*
- Textured, silky nipple for soothing comfort
- Designed for natural oral development

##### Sterilisation & carrying case

- Sterilizing in the microwave is safe and easy
- Simple steps for reliable sterilization in 3 minutes

##### Quality assured

- Made at our award-winning site in the UK

# PHILIPS

## Features

### Flexible, soft shield




Delicate skin needs extra care. Our FlexiFit technology lets this shield follow the natural curves of your baby's face for a comfortable fit. Your little one will experience fewer skin marks and less irritation.

### Rounded shield



Our rounded shield minimizes pressure for comfortable soothing that's gentle on your little one's cheeks.

### Extra airflow



Skin needs to breathe, especially your little one's. Our shield has 6 air holes for extra airflow, designed to reduce skin irritation.

### 98% accept silicone nipples



Babies know what they like! We asked moms how their little ones respond to our textured silicone nipples and 98% said that their baby accepts the Philips Avent ultra soft pacifiers.

### Silky textured nipple



This pacifier's silky, textured silicone nipple soothes and comforts your little one.

### Orthodontic nipple



Our collapsible silicone nipple has a symmetrical shape that respects your baby's palate, teeth and gums as he or she grows.

### Sterilizer / carrying case



With our carry case you can sterilize your baby's pacifier in the microwave and store hygienically when you're on the move.

### Case makes sterilizing easy



Place well-rinsed pacifiers into the case, add 25ml/9fl.oz. of water and microwave at 750-1000W for 3 minutes. Leave to cool for 5 minutes, pour water away and they're ready to go.

### Made in the UK



You can be sure that your little one's comfort is in good hands. This pacifier was made at our award-winning site in the UK.\*

## Specifications

Please note that this is a pre-sales leaflet. The contents of this leaflet reflect the best of our knowledge per date and country mentioned above. The contents of this leaflet are subject to change without notice. Philips does not accept any liability as to the contents of this leaflet.

### What is included

ultra soft pacifier 2 pcs

BPA-free Yes

### Accessories included

Sterilizer/carrying case Yes

### Hygiene

Easy to clean Yes  
Can be sterilized Yes  
Dishwasher safe Yes

### Safety

Safety ring handle Yes

- \* Developed with healthcare professionals and moms
- \* 2016-2017 US consumer tests show an average of 98% teat acceptance of the textured Philips Avent teat used in our ultra air and ultra soft soothers
- \* 2014 Manufacturer of the Year
- \* For hygiene reasons, replace pacifiers after 4 weeks of use
- \* No 1 global pacifier brand
- \* Our range supports moms and babies at every development stage
- \* 85% of moms surveyed perceive this pacifier feels softer than eight comparable models of other leading brands, independent research, USA, Feb 2017



data subject to change  
2019, July 20  
Version: 10.0.1  
UPC: 00075020065865

© 2019 Koninklijke Philips N.V.  
All rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

[www.philips.com](http://www.philips.com)