

PHILIPS

AVENT


ultra soft pacifier

ultra soft and flexible

6-18m

Orthodontic & BPA-Free

2-pack



SCF222/23



The softest soother for your baby's sensitive skin Flexible shield for a comfortable fit

Take care of sensitive skin with the Philips Avent ultra soft pacifier. Our super soft*, flexible shield follows the curves of your baby's cheeks, leaving fewer marks and less irritation for more comfortable soothing.

Unique soft and flexible shield

- Soft, flexible shield for fewer skin marks & less irritation
- Rounded shield for extra comfort every day

Comfortable soft & silky nipple

- Loved by babies worldwide*
- Textured, silky nipple for soothing comfort
- Designed for natural oral development

Sterilize and store in one handy case

- Sterilize and store in one handy case



asimpleswitch.com

Highlights

Flexible, soft shield



Sensitive skin needs extra care. Our shield technology lets this soother follow the natural curves of your baby's face for a comfortable fit. Your little one will experience fewer skin marks and less irritation.

Rounded shield



Our rounded shield minimizes pressure for comfortable soothing that's gentle on your little one's cheeks.

98% nipple acceptance*



When we asked moms how their little ones respond to our textured silicone nipples, on average 98% said that their baby accepts Philips Avent ultra soft and ultra air pacifiers.

Silky textured nipple



This pacifier's silky, textured silicone nipple soothes and comforts your little one.

Orthodontic nipple



Our collapsible silicone nipple has a symmetrical shape that respects your baby's palate, teeth and gums as he or she grows.

Easy sterilization & storage



The ultra soft travel case doubles as a sterilizer, so all you need to do is add some water and pop it in the microwave. Then rest easy that it's clean for the next use.

Specifications

What is included

ultra soft pacifier: 2 pcs

Accessories included

Sterilizer/carrying case

Safety

Safety ring handle
BPA-free

Hygiene

Easy to clean
Can be sterilized
Dishwasher safe



[asimpleswitch.com](https://www.asimpleswitch.com)

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.



* Developed with healthcare professionals and moms
* 2016-2017 US consumer tests show an average of 98% teat acceptance of the textured Philips Avent teat used in our ultra air and ultra soft soothers
* For hygiene reasons, replace pacifiers after 4 weeks of use
* No 1 global pacifier brand
* Our range supports moms and babies at every development stage
* 85% of moms surveyed perceive this pacifier feels softer than eight comparable models of other leading brands, independent research, USA, Feb 2017