

Philips
Nutrition starter pack



SCF230

Everything you need to prepare your baby's food,

all bottles and jars are sterile and safe now

Everything you need to quickly and confidently prepare your baby's food, with the added assurance that all bottles and jars are sterile and safe.

Baby's food never too hot

- Temperature control system
- The baby's food never gets too hot
- Graduated temperature control
- Safe cool wall due to double insulation

Fast & reliable sterilising

- Convenient steriliser for use in your microwave oven
- Kills all household bacteria in minutes using tap water
- Sterilises up to 6 regular bottles and 4 wide-neck bottles
- No chemicals needed

Quick & easy food preparation

- Bottle and jar lifter
- Click-step temperature setting
- Fits all bottle and jar types

PHILIPS
sense and simplicity

Specifications

Technical Specifications

- Heating time: 8-10 minute(s)
- Sterilising time: 8-12 minute(s)
- Safety Classification: Class 2

Convenience

- Automatic shut-off

- Fits all bottles: All baby bottles & jar types
- Power on indication
- Ready indicator

Power

- Bottle Warmer: 220/240 V

Highlights

Overheat protection

To avoid the danger of accidentally burning your baby's mouth, the automatic switch-off ensures that the food does not get too hot.

PTC heating element

To avoid the danger of accidentally burning your baby's mouth, the automatic switch-off ensures that the food does not get too hot.

Double insulation

With your safety and convenience in mind, the double insulated bottle warmer is designed to be used in every room in the house, not just where an earthed mains socket is available. There are no exposed metal parts that can become live.

Disinfects completely

Protect your newborn baby from any harmful household bacteria or germs by disinfecting her bottles, teats, soothers and small plastic toys.

Holds 4-6 bottles

Convenient use of the steriliser means having the disinfected bottles ready when you need them. This means fitting all bottles needed for a day's feeding into the steriliser in one go.

No chemicals needed

Chemical residues in baby's food are a nightmare to think about. Therefore, items are sterilised using very hot steam.

Fits all bottle and jar types

The product is designed in such a way that bottles and jars of all kinds of sizes and shapes fit in.



Issue date 2009-06-10

Version: 1.0

© 2009 Koninklijke Philips Electronics N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips
Electronics N.V. or their respective owners.

www.philips.com