

# PHILIPS

## AVENT

### ultra air pacifier

For baby's sensitive skin

6-18m

Orthodontic & BPA-Free

2-pack



SCF244/22



## A light, breathable pacifier for sensitive skin

Extra-large air holes let your baby's skin breathe

Let your little one's skin breathe with the Philips Avent ultra air pacifier. Featuring 4 extra-large air holes, the lightweight shield is designed for maximum air flow so your baby's sensitive skin stays drier while soothing.

#### **A light, breathable pacifier for sensitive skin**

- 4 extra large air holes
- Skin stays drier while soothing
- Rounded edges for comfortable soothing

#### **Comfortable soft & silky nipple**

- 98% of babies accept our ultra air pacifier\*
- Textured, silky nipple for soothing comfort
- Designed for natural oral development

#### **Sterilisation & carrying case**

- Sterilizing in the microwave is safe and easy
- Simple steps for reliable sterilization in 3 minutes



[asimpleswitch.com](http://asimpleswitch.com)

# Highlights

## 4 extra large air holes



ultra air is designed to allow maximum air flow, so your baby's sensitive skin can breathe.

## Skin stays drier



Your little one's sensitive skin stays drier while soothing, thanks to this pacifier's breathable design with maximum air flow.

## Comfortable, rounded edges



The ultra air shield is lightweight and features rounded edges for your baby's comfort.

## Nipple acceptance



Babies know what they like! We asked moms how their little ones respond to our textured silicone nipples and on average 98% said that their baby accepts the Philips Avent ultra soft and ultra air pacifiers.

## Silky textured nipple



This pacifier's silky, textured silicone nipple soothes and comforts your little one.

## Orthodontic nipple



Our collapsible silicone nipple has a symmetrical shape that respects your baby's palate, teeth and gums as he or she grows.

## Sterilizer / carrying case



With our carry case you can sterilize your baby's pacifier in the microwave and store hygienically when you're on the move.

## Case makes sterilizing easy



Place well-rinsed pacifiers into the case, add 25ml/1 fl. oz. of water and microwave at 750-1000W for 3 minutes. Leave to cool for 5 minutes, drain water and they're ready to go.



[asimpleswitch.com](http://asimpleswitch.com)

## Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

# Specifications

## What is included

ultra air pacifier: 2 pcs

## Accessories included

Sterilizer/carrying case

## Safety

Safety ring handle

BPA-free

## Hygiene

Easy to clean

Can be sterilized

Dishwasher safe

