

PHILIPS

AVENT

ultra air pacifier 0-6m, 2 pack

For baby's sensitive skin

0-6m

Orthodontic & BPA-Free

2-pack



SCF245/20



A light, breathable pacifier for sensitive skin

Extra-large air holes let your baby's skin breathe

Let your little one's skin breathe with the Philips Avent ultra air pacifier. Featuring 4 extra-large air holes, the lightweight shield is designed for maximum air flow so your baby's sensitive skin stays drier while soothing.

A light, breathable pacifier for sensitive skin

- 4 extra large air holes
- Skin stays drier while soothing
- Rounded edges for comfortable soothing

Comfortable soft & silky nipple

- 98% of babies accept our ultra air pacifier*
- Textured, silky nipple for soothing comfort
- Designed for natural oral development

Sterilisation & carrying case


- Sterilizing in the microwave is safe and easy
- Simple steps for reliable sterilization in 3 minutes



asimpleswitch.com

Highlights

4 extra large air holes



ultra air is designed to allow maximum air flow, so your baby's sensitive skin can breathe.

Skin stays drier



Your little one's sensitive skin stays drier while soothing, thanks to this pacifier's breathable design with maximum air flow.

Comfortable, rounded edges



The ultra air shield is lightweight and features rounded edges for your baby's comfort.

Nipple acceptance



Babies know what they like! We asked moms how their little ones respond to our textured silicone nipples and on average 98% said that their baby accepts the Philips Avent ultra soft and ultra air pacifiers.

Silky textured nipple



This pacifier's silky, textured silicone nipple soothes and comforts your little one.

Orthodontic nipple



Our collapsible silicone nipple has a symmetrical shape that respects your baby's palate, teeth and gums as he or she grows.

Sterilizer / carrying case



With our carry case you can sterilize your baby's pacifier in the microwave and store hygienically when you're on the move.

Case makes sterilizing easy



Place well-rinsed pacifiers into the case, add 25ml/1 fl. oz. of water and microwave at 750-1000W for 3 minutes. Leave to cool for 5 minutes, drain water and they're ready to go.



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

Specifications

What is included

ultra air pacifier: 2 pcs

Accessories included

Sterilizer/carrying case

Safety

Safety ring handle
BPA-free

Hygiene

Easy to clean
Can be sterilized
Dishwasher safe

