

# PHILIPS

## AVENT

Bottle to Cup  
Trainer Kit

Natural

150 ml/5 oz

4m+

3 hole Medium Flow



SCF251/00



## Ease your baby's transition to a drinking cup

Our new Natural bottle-to-cup trainer kit helps your baby transition from bottle drinking to his first cup easily. With the soft-touch grips, your baby will be able to hold the bottle independently while drinking from the familiar teat.

### Ideal first step to a drinking cup

- Soft touch trainer handles for little hands
- Advanced anti-colic system with innovative twin valve
- Unique Petals for a soft, flexible teat without collapse

### Convenient for you and your baby

- Interchangeable nipple and soft spout
- All parts are dishwasher safe for convenience

### Other benefits

- Compatible with Philips Avent Natural and Classic+ bottles
- This cup is made from BPA-free material



[asimpleswitch.com](http://asimpleswitch.com)

# Highlights

## Unique comfort petals



Petals inside the teat increase softness and flexibility without teat collapse. Your baby will enjoy a more comfortable and contented feed.

## Advanced anti-colic valve



Innovative twin valve designed to reduce colic and discomfort by venting air into the bottle and not baby's tummy.

## Cup is dishwasher safe



All parts are dishwasher safe for convenience

## BPA-Free



This Philips Avent cup is made from BPA-free material.

## Compatibility across the range



The new Philips Avent cup is compatible with the Philips Avent Natural and Classic+ bottle range (excluding the glass bottle) and all other Philips Avent cups (excluding the Natural Drinking Cup).

## Ergonomically designed trainer



Trainer handles help your toddler hold the cup and drink independently. These handles are not only shaped for little hands to grab easily, they are also rubberised for a non-slip grip.

## Interchangeable



The sip, no-drip soft spouts are interchangeable with natural nipples. This gives you flexibility to move from bottle to cup whenever your child is ready.



[asimpleswitch.com](http://asimpleswitch.com)

## Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas — Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

# Specifications

## What is included

Bottle to 1st trainer cup: 1 pcs  
Snap on hygienic cap: 1 pcs  
Soft, Medium Flow Teat: 1 pcs  
Size: 150 ml/5 oz

## Development stages

Stages: 4 months +

## Country of origin

Indonesia

## Material

Bottle: BPA free\*, Polypropylene  
Teat: BPA free\*, Silicone

## Weight and dimensions

No. of F-boxes in A-box: 6 pcs

