

PHILIPS

AVENT

Bottle to Cup
Trainer Kit

Natural

4m+



SCF251/04



Ease your baby's transition to a drinking cup

This Natural trainer kit eases your baby into their first cup. The comfort grips helps your baby hold the bottle independently while drinking from the familiar nipple. When ready, switch to the silicone spout to complete the transition.

Ideal first step to a drinking cup

- Soft touch trainer handles for little hands
- Advanced anti-colic system with innovative twin valve
- Unique Petals for a soft, flexible nipple without collapse

Convenient for you and your baby

- Compatible with Philips Avent bottles and cups
- All parts are dishwasher safe for convenience

Other benefits

- This cup is made from BPA free material
- Compatible with Philips Avent bottles and cups



asimpleswitch.com


Highlights

Ergonomic trainer handles



Trainer handles help your toddler to hold the cup and drink independently. These handles are shaped for little hands to grab easily and they are also rubberized for a non-slip grip.

Unique comfort petals



Petals inside the nipple increase softness and flexibility without nipple collapse. Your baby will enjoy a more comfortable and contented feed.

Advanced anti-colic valve



Innovative twin valve designed to reduce colic and discomfort by venting air into the bottle and not baby's tummy.

Cup is dishwasher safe



All parts are dishwasher safe for convenience

BPA Free



This Philips Avent cup is made from BPA free material.

Philips Avent compatibility



All Philips Avent bottles and cups are compatible excluding the glass bottles and the Grown-up cups/My First Big Kid Cups. So you can mix and match to create the perfect cup, suiting your toddler's individual development needs.

Philips Avent compatibility



All Philips Avent bottles and cups are compatible excluding the glass bottles and the Grown-up cups/My First Big Kid Cups. So you can mix and match to create the perfect cup, suiting your toddler's individual development needs.



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

Specifications

What is included

Bottle to 1st trainer cup: 1 pcs
Snap-on hygienic cap: 1 pcs
Soft, Medium Flow Nipple: 1 pcs
Silicone spout: 1 pcs

Development stages

Stage: 4 months +

Country of origin

Indonesia

Material

Bottle: BPA free*, Polypropylene
Nipple: BPA free*, Silicone

Weight and dimensions

F-box dimensions: 198 x 125 x 9 mm
No. of F-boxes in A-box: 6 pcs
Product weight: 0.11 kg

Accessories

Replacement part: Replacement spout SCF252

