

PHILIPS

AVENT

Microwave Steriliser
Starter Set

SCF271/02



Ultra fast and convenient

Sterilises 6 bottles in 2 minutes

Philips Avent Microwave Steam Steriliser's lightweight, compact design makes it ideal for use in and out of the home. Contents remain sterile for up to 24 hours if lid is unopened.

Benefits and features:

- Ultra fast, easy to use
- Compact and lightweight
- Ideal for home and travel
- Large capacity
- Safe, easy handling
- Effective sterilisation

Highlights

Ultra fast, easy to use



Just add water, load and place in the microwave for as little as 2 mins. 2 mins at 1100-1850 Watt, 4 mins at 850-1000 Watt, 6 mins at 500-850 Watt.

Compact and lightweight



Fits most microwaves. Convenient for travel.

Large capacity



Holds up to six 260-ml/9-oz Avent baby bottles or two Philips Avent Breast Pumps and two Avent baby bottles.

Safe, easy handling



Clips close lid securely and side grips stay cooler to aid safe handling.

Effective sterilisation



Based on the hospital principle, the intense heat of the steam destroys harmful bacteria.

Specifications

Weight and dimensions

Weight: 740 g

Dimensions: 166 (H), 280 (W), 280 (L) mm

Country of origin

England

What is included

Microwave steam steriliser: 1 pcs

Airflex Feeding Bottle (260 ml/9 oz): 2 pcs

Extra-soft slow flow teat: 2 pcs

Tongs: 1 pcs

Development stages

Stages: 0-6 months

