

PHILIPS

AVENT

Microwave Steam
Steriliser

SCF271/06



Ultra fast and convenient

Sterilises 6 bottles in 2 minutes*

Philips Avent Microwave Steam Steriliser's lightweight, compact design makes it ideal for use in and out of the home. Contents remain sterile for up to 24 hours if lid is unopened.

Ultra fast, easy to use

- Just add water, load and place in the microwave

Compact and lightweight

- Convenient for travel. Fits most microwaves.

Safe, easy handling

- Contents remain sterile for up to 24 hours if unopened

Effective sterilisation

- Extra protection through natural steam sterilisation

Safe easy handling

- Side grips close the lid securely

Specifications

Weight and dimensions

Dimensions: 166 (H), 280 (W), 280 (L) mm

Weight: 740 g

Country of origin

England

What is included

Microwave steam steriliser: 1 pcs

Feeding Bottle (260 ml/9 oz): 1 pcs

Feeding Bottle (125 ml/4 oz): 1 pcs

Extra soft Slow Flow teat: 1 pcs

Extra soft Newborn Flow teat: 1 pcs

Newborn Soother: 1 pcs

Tongs: 1 pcs

Development stages

Stages: 0–6 months

