

PHILIPS

AVENT

Electric Steam
Sterilizer

110V

SCF274/30



Fast and easy to use

One touch operation

This steriliser is simple and straightforward to use right from the start. At the touch of a button it will sterilise your bottles and accessories in only 8 minutes.

Large capacity

- Holds up to six Philips Avent Bottles

Easy to use

- Just add water, load and switch on

Fast

- Contents are sterile and ready to use in approx. 8 minutes

Safe, easy handling

- Contents remains sterile for up to 6 hours if unopened

Effective sterilization

- Extra protection through natural steam sterilization

Highlights

Holds up to 6 bottles



The sterilizer has a smart design that takes up little space in your kitchen, yet still fits six Philips Avent Bottles or two Philips Avent Breast Pumps. The two baskets inside can clip together to form a dishwasher basket, which makes pre-cleaning of smaller items like soothers and nipples much easier.

Just add water



The Philips Avent Electric Sterilizer is very easy to use. Just add water, load the sterilizer with the contents you would like to sterilize and switch it on.

Contents sterile in 8 minutes



With the Philips Avent Electric Sterilizer you can sterilize contents in as little as 8 minutes

Sterile contents for 6 hours



Once the sterilization cycle has been completed, the contents inside the Philips Avent Electric Sterilizer will remain sterile for up to 6 hours if the lid is unopened.

Natural steam sterilization

Sterilizing is all about protecting your baby from particularly harmful milk bacteria until his immune system is strong enough. The Philips Avent Sterilizer uses the hospital method of steam sterilization, which is quick, easy, and efficient, with no chemicals involved.

Specifications

Weight and dimensions

Dimensions: 313 (H), 235 (W), 235 (L) mm

Weight: 1.575 kg

Country of origin

England

What is included

Tongs: 1 pcs

Power

Voltage: 110 V

Development stages

Stage: 0 - 6 months

