

PHILIPS

AVENT

Electric Steam
Sterilizer

220-240V

SCF274/31



Sterilizes 6 bottles in 8 minutes

Fast and easy to use

Sterile contents for up to 6 hours

Fast

- Contents are sterile and ready to use in approx. 8 minutes

Safe, easy handling

- Contents remains sterile for up to 6 hours if unopened

Easy to use

- Just add water, load and switch on

Effective sterilization

- Extra protection through natural steam sterilization

Large capacity

- Holds up to six Philips AVENT Bottles

Specifications

Weight and dimensions

Dimensions: 313 (H), 235 (W), 235 (L) mm

Weight: 1.575 kg

Country of origin

England

What is included

Tongs: 1 pcs

Power

Voltage: 220 - 240 V

Development stages

Stages: 0 - 6 months

