

PHILIPS

AVENT

24-Hour Steam Steriliser

SCF276/01



Keeps contents sterile night and day

Add or remove items at any time

The easiest way to sterilise is with the Philips Avent 24 Hour Steriliser. It keeps contents sterile for 24 hours and the digital display tells you when it is time to re-sterilise. Includes 2 Avent Feeding Bottles and handy accessories.

Automatic 24-hour sterilisation

- Sterile contents, even if you open the lid

Easy to use

- Just add water once a day

Keeps you informed

- Advanced digital display and sound alerts

Fast

- Contents are sterile and ready to use in approx. 6 min.

Large capacity

- Holds up to six Philips Avent Bottles


Highlights

Keeps contents sterile



This advanced steriliser makes sure you always have a sterile bottle on hand by keeping contents sterile for 24 hours. It will automatically re-sterilise the contents every 6 hours. With the handy pause button you can also add or remove items at any time and the 24-Hour Steriliser will automatically re-sterilise the contents.

Advanced digital display



The display and sound alerts tell you when sterilisation is complete and how long items will remain sterile. The display also lets you know when you have filled the steriliser with too little or too much water, so you can always be sure that sterilising has been completed accurately and reliably.

Just add water once a day



Aside from its advanced functionality, the 24-Hour Steriliser is simple and straightforward to use. Just add water, load the steriliser with

your bottles, breast pump or accessories and select the mode. Mode 1 lets you run a single sterilising cycle, which is useful when you want to take all the items out in one go. Mode 2 is for 24-hour sterilising, which is very convenient when you want to leave the items in the steriliser overnight or when you want to take out only one bottle for every feed. No matter what your routine will be, this steriliser will make the daily sterilising chore easier.

Holds up to six bottles



The steriliser has a smart design that takes up little space in your kitchen, yet still fits six Philips Avent Bottles or two Philips Avent Breast Pumps. The two baskets inside can clip together to form a dishwasher basket, which makes pre-cleaning of smaller items like soothers and teats much easier.

Sterile contents in 6 minutes



The advantage of steam sterilisation is that it is much quicker (and safer) than the traditional method of boiling bottles in a pan on the stove. The 24-Hour Steriliser takes as little as 6 minutes to sterilise a full load of 6 bottles, teats and dome caps.

Natural steam sterilisation



Sterilising is all about protecting your baby from particularly harmful milk bacteria until his immune system is strong enough. The Philips Avent Steriliser uses the hospital method of steam sterilisation, which is quick, easy and efficient, with no chemicals involved.



Practical Parenting Award

Every year Practical Parenting magazine rates baby products in every baby category. Finalist in the 2009 Practical Parenting Awards



Prima Baby & Pregnancy Readers' Awards

Prima Baby asks its readers to vote on the best baby products on the market. Voted 'Best Buy' by Prima Baby readers

Specifications

Weight and dimensions

Dimensions: 320 (H), 230 (W), 230 (L) mm

Weight: 1.47 kg

Country of origin

England

What is included

Tongs: 1 pcs

Newborn Soother: 1 pcs

Measuring jug: 1 pcs

Bottle and teat brush: 1 pcs

Extra-soft newborn flow teat: 2 pcs

Extra Durable Feeding Bottle: 125 ml/4 oz,

260 ml/9 oz

Power

Voltage: 220-240 V

Development stages

Stages: 0-6 months

