

# PHILIPS

## AVENT

### 3-in-1 Electric Steam Steriliser

SCF284/02



## The easiest way to sterilise

### Adjustable size saves kitchen space

Our new electric steam steriliser has been designed to make sterilising as simple as possible. With its adjustable size, it takes up less kitchen space and the baskets fit both wide and narrow neck bottles.

#### **3-in-1 steriliser**

- Adjustable size

#### **Fits all bottle types**

- Suitable for wide neck and narrow neck bottles

#### **Extra speed and safety**

- 6-minute cycle with automatic shut-off

#### **Effective sterilisation**

- Natural steam sterilisation kills 99.9% of harmful germs

#### **BPA-Free**

- This steriliser is made from BPA-Free materials

# Highlights

## 3-in-1 steriliser



Due to its adjustable size, this 3-in-1 steriliser perfectly fits the items you would like to sterilise and always takes up the least amount of kitchen space. The 3-in-1 functionality gives you: 1) Small size configuration for sterilising soothers 2) Medium size configuration for sterilising breast pumps, toddler plates, knives and forks 3) Large size configuration for sterilising 6 bottles

## Fits all bottle types



The steriliser has been designed to fit both wide neck and narrow neck bottles as well as different breast pumps and accessories. This gives you the freedom to choose the bottles and accessories that best suit you and your baby. At full capacity the steriliser holds six 11 oz/330 ml Philips AVENT bottles.

## 6-minute cycle and auto shut-off



Your bottles and accessories will be sterilised in just 6 minutes. Once sterilisation is complete, the device will shut off automatically, ensuring extra safety and less energy consumption.

## Kills 99.9% of harmful germs



Sterilising protects your baby from particularly harmful bacteria until his immune system is strong enough. Philips AVENT Sterilisers use the hospital method of steam sterilisation, which is quick, easy and effective, with no chemicals involved. Contents stay sterile for 24 hours if the lid is unopened.

## BPA-Free



This steriliser is made from BPA-Free materials

# Specifications

## Country of origin

Turkey

## What is included

Tongs: 1 pcs

## Power

Voltage: 220 V

## Development stages

Stages: 0–6 months

