

PHILIPS

AVENT

3-in-1 Electric
Steam Sterilizer

Electric Sterilizer

SCF284/03



The easiest way to sterilize

Adjustable size saves kitchen space

Our new electric steam sterilizer has been designed to make sterilizing as simple as possible. With adjustable size, it takes up less kitchen space and the baskets fit both wide and narrow neck bottles.

3-in-1 sterilizer

- Adjustable size

Fits all bottle types

- Suitable for wide neck and narrow neck bottles

Extra speed and safety

- 6 minute cycle with automatic shut-off

Effective sterilization

- Natural steam sterilization kills 99.9% of harmful germs

BPA Free

- This sterilizer is made from BPA Free materials



Highlights

3-in-1 sterilizer



Due to its adjustable size, this 3-in-1 Philips AVENT sterilizer perfectly fits the items you would like to sterilize and always takes up the least amount of kitchen space. The 3-in-1 functionality gives you: 1) Small size configuration for sterilizing soothers 2) Medium size configuration for sterilizing breast pumps, toddler plates, knives, and forks 3) Large size configuration for sterilizing 6 baby bottles

Fits all bottle types



The sterilizer has been designed to fit both wide neck and narrow neck bottles as well as different breast pumps and accessories. This gives you the freedom to choose the bottles and accessories that best suit you and your baby. At full capacity the sterilizer holds six 11oz / 330ml Philips AVENT bottles.

6 minute cycle & auto shut-off



Your bottles and accessories will be sterilized in just 6 minutes. Once sterilization is complete the device will shut off automatically, ensuring extra safety and less energy consumption.

Kills 99.9% of harmful germs



Sterilizing protects your baby from particularly harmful bacteria until his immune system is strong enough. Philips AVENT Sterilizers use the hospital method of steam sterilization, which is quick, easy, and effective, with no chemicals involved. Contents stay sterile for 24 hours if the lid is unopened.

BPA Free



This sterilizer is made from BPA Free materials

Specifications

Country of origin

Turkey

What is included

Tongs: 1 pcs

Power

Voltage: 220 V

Development stages

Stages: 0 - 6 months

Compatibility

Philips-AVENT range compatible

Material

Polypropylene

Weight and dimensions

Dimensions: 290 x 160 x 350 (w x d x h) mm

Weight: 1.5 kg



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.