

# PHILIPS

## AVENT

### 3-in-1 electric steam steriliser

Kills 99.9% of harmful germs

Sterilises in 6 minutes

Fits 6 Philips Avent bottles

Adjustable 3-in-1 design

SCF285/01



## Convenient and effective sterilisation

### Flexible, easy loading

With its adjustable size, the Philips Avent 3-in-1 electric steam steriliser takes up the least amount of kitchen space, yet perfectly fits the items you want to sterilise, whether it's a few small items or a full load.

#### Space-saving 3-in-1 design

- 3-in-1 modular steriliser design

#### Effective sterilising using natural steam

- Natural steam sterilising kills 99.9% of harmful germs
- Sterile for up to 24 hours if the lid is unopened
- Fast cycle and auto shut off function

#### Fits various bottle types and other products

- Spacious design
- Sterilises various bottles, breast pumps and accessories

#### Easy-to-clean and safe sterilising

- Easy-to-clean and safe sterilising

## Highlights

### Flexible 3-in-1 steriliser



The steriliser's unique modular design makes it quick and easy to fit and organise your bottles and accessories. It also takes up minimal kitchen space.

### Natural steam sterilising



The steriliser uses natural steam to sterilise baby bottles and other products, killing 99.9% of harmful germs without using chemicals. Perfect for giving you peace of mind that all your baby bottles and other products are sterile.

### Sterile for up to 24 hours



The steriliser will keep its contents—baby bottles, breast pumps, etc.—sterile for up to 24 hours if the lid is unopened.

### Fast sterilisation



The sterilising cycle lasts just 6 minutes. After this time, the steriliser automatically shuts off, giving you extra speed and safety.

### Spacious design



Holds up to six 330 ml / 11 oz Philips Avent Classic and Natural bottles to be sterilised at once.

### Open design



With an open design, it's easy to keep the heating plate clean, so you can sterilise with clean steam every time.

### Works across bottle types

The steriliser is suitable for use with both standard-neck and wide-neck baby bottles. It also fits other baby products like breast pumps and accessories.

## Specifications

### Weight and dimensions

Dimensions: 290 x 160 x 150 (w x d x h) mm  
Weight: 1.5 kg

### Country of origin

Made in: Turkey

### What is included

Electric steam steriliser: 1 pcs  
Classic Bottle 4 oz: 1 pcs  
Newborn soother: 1 pcs  
Tongs: 1 pcs

### Power

Voltage: 220–240 V

### Development stages

Stages: 0–6 months

### Compatibility

Philips Avent range-compatible

### Material

Polypropylene

### Technical specifications

Power consumption: 650 W

Voltage: 50–60Hz

Safety Classification: Class 1

