

PHILIPS

AVENT

3-in-1 Electric Steam Steriliser

Electric Steriliser

Sterilizes in 6 minutes

Fits 6 Philips Avent bottles

Adjustable 3-in-1 design

SCF285/02



The easiest way to sterilise

Adjustable steriliser, saves kitchen space

Our new Philips Avent electric steam steriliser has been designed to make sterilising as simple as possible. With adjustable size, it takes up less kitchen space and the baskets fit both wide and narrow neck bottles.

Space-saving 3-in-1 design

- 3-in-1 modular sterilizer design

Effective sterilizing using natural steam

- Natural steam sterilizing kills 99.9% of harmful germs
- Sterile for up to 24 hours if the lid is unopened
- Fast cycle and auto shut off function

Fits various bottle types and other products

- Spacious design
- Sterilizes various bottles, breast pumps & accessories

Easy to clean and safe sterilizing

- Easy to clean and safe sterilizing

Highlights

Flexible 3-in-1 sterilizer



The sterilizer's unique modular design enables you to fit the bottles and accessories flexibly, and organize them easily. Thus loading and unloading are very convenient. It also takes up minimal kitchen space.

Natural steam sterilizing



The sterilizer uses natural steam to sterilize baby bottles and other products, killing 99.9% of harmful germs without using chemicals. Perfect for giving you peace of mind that all your baby bottles and other products are sterile.

Sterile for up to 24 hours



The sterilizer will keep its contents – baby bottles, breast pumps etc – sterile for up to 24 hours if the lid is unopened.

Fast sterilization



The sterilizing cycle lasts just 6 minutes. After this time, the sterilizer automatically shuts off, giving you extra speed and safety.

Spacious design



Holds up to six 330 ml / 11oz Philips Avent Classic and Natural bottles to be sterilized at once.

Open design



With our open design you can easily keep the heating plate clean, allowing you to sterilize with clean steam every time.

Works across bottle types

The sterilizer is suitable for use with both standard neck and wide neck baby bottles. It also fits other baby products like breast pumps and accessories.

Specifications

Weight and dimensions

Dimensions: 290 x 160 x 150 (w x d x h) mm
Weight: 1.5 kg

Country of origin

Made in: Turkey

What is included

Electric steam sterilizer: 1 pcs
Classic Bottle 4oz: 1 pcs
Newborn pacifier: 1 pcs
Tongs: 1 pcs

Power

Voltage: 220–240 V

Development stages

Stages: 0 – 6 months

Compatibility

Philips-Avent range compatible

Material

Polypropylene

Technical specifications

Power consumption: 650 W

Voltage: 50–60Hz

Safety Classification: Class 1

