

# PHILIPS

## AVENT

4-in-1 Electric  
Steam Sterilizer

Electric Sterilizer

SCF286/02



## The easiest way to wash and sterilize

### Adjustable size saves kitchen space

The AVENT electric steam sterilizer SCF286/02 has been designed to make sterilizing as simple as possible. The 4-in-1 sterilizer has an integrated dishwasher basket for hassle free pre-cleaning and an advanced display to keep you informed.

#### **4-in-1 sterilizer**

- Adjustable size and integrated dishwasher basket

#### **Keeps you informed**

- Advanced display with progress and cool-down indication

#### **Fits all bottle types**

- Suitable for wide neck and narrow neck bottles

#### **Extra speed and safety**

- 6 minute cycle with automatic shut-off

#### **Effective sterilization**

- Natural steam sterilization kills 99.9% of harmful germs

# Highlights

## 4-in-1 sterilizer



Due to its adjustable size, this 4-in-1 Philips AVENT sterilizer perfectly fits the items you would like to sterilize and always takes up the least amount of kitchen space. The 4-in-1 functionality gives you: 1) Small size configuration for sterilizing soothers 2) Medium size configuration for sterilizing several bottles, toddler plates, knives, and forks 3) Large size configuration for sterilizing a full load of up to 6 bottles 4) Integrated dishwasher basket for hassle free pre-cleaning of baby bottles and accessories

## Advanced display



The advanced display keeps you informed throughout the sterilization cycle and lets you know when the bottles are cooled down.

## Fits all bottle types



The sterilizer has been designed to fit both wide neck and narrow neck bottles as well as different breast pumps and accessories. This

gives you the freedom to choose the bottles and accessories that best suit you and your baby. At full capacity the sterilizer holds six 11oz / 330ml Philips AVENT bottles.

## 6 minute cycle & auto shut-off



Your bottles and accessories will be sterilized in just 6 minutes. Once sterilization is complete the device will shut off automatically, ensuring extra safety and less energy consumption.

## Kills 99.9% of harmful germs



Sterilizing protects your baby from particularly harmful bacteria until his immune system is strong enough. Philips AVENT Sterilizers use the hospital method of steam sterilization, which is quick, easy, and effective, with no chemicals involved. Contents stay sterile for 24 hours if the lid is unopened.

## BPA Free



This sterilizer is made from BPA Free materials

# Specifications

## Country of origin

Turkey

## What is included

Tongs: 1 pcs

## Power

Voltage: 220 V

## Development stages

Stages: 0 - 6 months

## Compatibility

Philips-AVENT range compatible

## Material

Polypropylene

## Weight and dimensions

Dimensions: 290 x 160 x 350 (w x d x h) mm

Weight: 1.5 kg

