

# PHILIPS

## AVENT

### 4-in-1 Electric Steam Steriliser

Electric Steriliser

Sterilises in 6 minutes

Fits 6 Philips Avent bottles

Adjustable 4-in-1 design

SCF287



## The easiest way to sterilise

### Adjustable size and dishwasher basket

Our new electric steam steriliser (also spelt sterilizer) has been designed to make sterilising bottles as simple as possible. The 4-in-1 steriliser has an integrated dishwasher basket for hassle free pre-cleaning and an advanced display to keep you informed.

#### Space-saving 4-in-1 design

- 4-in-1 modular sterilizer design

#### Fast, effective sterilisation using natural steam

- Natural steam sterilizing kills 99.9% of harmful germs
- Sterile for up to 24 hours if the lid is unopened
- Fast cycle and auto shut off function
- Integrated display shows progress

#### Easy to clean and safe sterilising

- Easy to clean and safe sterilizing

#### Fits various bottle types and other products

- Spacious design
- Sterilises various bottles, breast pumps & accessories

# Highlights

## Flexible 4-in-1 sterilizer



The sterilizer's unique modular design enables you to fit the bottles and accessories flexibly, and organize them easily. Thus loading and unloading are very convenient. It also takes up minimal kitchen space. The included dishwasher basket keeps small items together in the dishwasher and allows you to transfer them all at once into the sterilizer.

## Natural steam sterilizing



The sterilizer uses natural steam to sterilize baby bottles and other products, killing 99.9% of harmful germs without using chemicals. Perfect for giving you peace of mind that all your baby bottles and other products are sterile.

## Sterile for up to 24 hours



The sterilizer will keep its contents - baby bottles, breast pumps etc - sterile for up to 24 hours if the lid is unopened.

## Fast sterilization



The sterilizing cycle lasts just 6 minutes. After this time, the sterilizer automatically shuts off, giving you extra speed and safety.

## Integrated display



The steriliser features an integrated display to show how the sterilising cycle is progressing. It also lets you know when the bottles have cooled down and are safe to remove.

## Open design



With our open design you can easily keep the heating plate clean, allowing you to sterilize with clean steam every time.

## Spacious design



Holds up to six 330 ml / 11oz Philips Avent Classic and Natural bottles to be sterilized at once.

## Works across bottle types

The steriliser is suitable for use with both standard neck and wide neck baby bottles. It also fits other baby products like breast pumps and accessories.

# Specifications

## Weight and dimensions

Dimensions: 290 x 160 x 150 (w x d x h) mm

Weight: 1.5 kg

## Country of origin

Made in: Turkey

## What is included

Electric steam sterilizer: 1 pcs

Classic Bottle 4oz: 1 pcs

Newborn pacifier: 1 pcs

Tongs: 1 pcs

## Power

Voltage: 220-240 V

## Development stages

Stages: 0 - 6 months

## Compatibility

Philips-Avent range compatible

## Material

Polypropylene

## Technical specifications

Power consumption: 650 W

Voltage: 50-60Hz

Safety Classification: Class 1

