

PHILIPS

AVENT

4-in-1 Electric
Steam Steriliser

Electric Steriliser

SCF287/02



The easiest way to sterilise

Adjustable size saves kitchen space

Our new electric steam steriliser has been designed to make sterilising bottles as simple as possible. The 4-in-1 steriliser has an integrated dishwasher basket for hassle-free pre-cleaning and an advanced display to keep you informed.

4-in-1 steriliser

- Adjustable size and integrated dishwasher basket

Keeps you informed

- Advanced display with progress and cool-down indication

Fits all bottle types

- Suitable for wide neck and narrow neck bottles

Extra speed and safety

- 6-minute cycle with automatic shut-off

Effective sterilisation

- Natural steam sterilisation kills 99.9% of harmful germs

Highlights

4-in-1 steriliser



Due to its adjustable size, this 4-in-1 Philips AVENT steriliser perfectly fits the items you would like to sterilise and always takes up the smallest amount of kitchen space. The 4-in-1 functionality gives you: 1) Small size configuration for sterilising soothers 2) Medium size configuration for sterilising several bottles, toddler plates, knives and forks 3) Large size configuration for sterilising a full load of up to 6 bottles 4) Integrated dishwasher basket for hassle-free pre-cleaning of baby bottles and accessories

Advanced display



The advanced display keeps you informed throughout the sterilisation cycle and lets you know when the bottles are cooled down.

Fits all bottle types



The steriliser has been designed to fit both wide neck and narrow neck bottles as well as different breast pumps and accessories. This gives you the freedom to choose the bottles

and accessories that best suit you and your baby. At full capacity the steriliser holds six 11 oz/330 ml Philips AVENT bottles.

6-minute cycle and auto shut-off



Your bottles and accessories will be sterilised in just 6 minutes. Once sterilisation is complete, the device will shut off automatically, ensuring extra safety and less energy consumption.

Kills 99.9% of harmful germs



Sterilising protects your baby from particularly harmful bacteria until his immune system is strong enough. Philips AVENT Sterilisers use the hospital method of steam sterilisation, which is quick, easy and effective, with no chemicals involved. Contents stay sterile for 24 hours if the lid is unopened.

BPA-Free



This steriliser is made from BPA-Free materials

Specifications

Weight and dimensions

Dimensions: 290 x 160 x 150 (w x d x h) mm

Weight: 1.5 kg

Country of origin

Turkey

What is included

Tongs: 1 pcs

Newborn Soother: 1 pcs

Classic Bottle 4 oz: 1 pcs

Power

Voltage: 220 V

Development stages

Stages: 0–6 months

Compatibility

Philips AVENT range compatible

Material

Polypropylene

