



Philips AVENT
Manual breast pump

Out and about set



SCF290/13

Express more. Quickly.

Enhanced efficiency from a manual breast pump

A manual breast pump offers you total control of the suction level and pumping rhythm. The unique Philips AVENT Manual Breast Pump assures you enhanced efficiency, being clinically proven to express more milk than a hospital electric pump.

Clinically proven effective

- Patented Let-down Massage Cushion

Flexible and portable

- Easy to assemble and use

Stimulate and express

- Soft massage cushions trigger natural let-down, like baby

Benefits and features:

- With insulated travel bag and cool packs

PHILIPS
AVENT

sense and simplicity

Specifications

What is included

- Manual breast pump: 1 pcs
- Breast milk container (125 ml/4 oz): 2 pcs
- Breast milk container (260 ml/9 oz): 2 pcs
- Extra soft, Newborn Flow Teat: 2 pcs
- Cool pack: 2 pcs
- Bottle stand/funnel cover: 2 pcs

- Microfibre insulated travel bag: 1 pcs
- Spare parts: 1 pcs

Country of origin

- England

Development stages

- Stages: 0–6 months

Highlights

With insulated travel bag and cool packs



Everything you need to maintain your milk supply and to express and store your breast milk when away from your baby.

Patented Let-down Massage Cushion



Five, soft petal massagers gently flex in and out as you pump and work with the vacuum to imitate your baby's suckling.

Soft massage cushions



Patented massage cushions flex in and out, replicating baby's suckling action and stimulating a fast, natural let-down



Issue date 2012-07-26

Version: 7.3.1

12 NC: 0000 000 00000

© 2012 Koninklijke Philips Electronics N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips
Electronics N.V. or their respective owners.

www.philips.com